

# Semi and fully automated immunostaining sample preparation platforms improve live leukocyte recovery, reproducibility, and flow cytometry data quality

## ABSTRACT


**BACKGROUND:** Limited innovation in automated cell and organelle sample preparation methodology limits the effectiveness of modern analytical methods, such as single-cell 'omics, flow and mass cytometry. These techniques traditionally rely on manual centrifugation-based protocols for cell washing and suspension preparation, hampering researchers' access to the reproducibility and scalability benefits of automation.

**METHODS:** We have developed a suite of cell suspension preparation systems that enable semi and full automation of cell washing protocols. These Laminar Wash™ technologies robustly, gently, and efficiently remove debris, dead cells, and unbound reagent using laminar flow and liquid handling robotics, rather than turbulent and harsh pelleting-plus-pipetting methods. Murine and humanized mouse peripheral blood mononuclear cells (PBMCs) and tumor infiltrating lymphocytes (TILs) were prepared and immunostained for flow cytometry analysis. Workflow improvements were assessed, as well as data quality by flow cytometry gating strategies isolating live cells and various lymphocyte subpopulations.


**RESULTS:** Adaptation of standard protocols to Laminar Wash automation typically improves repetitive immunostaining processes and workflows, in terms of reduced hands-on time and inter- and intra-operator variability. We demonstrate the superior live cell retention and reproducibility of Laminar Wash over centrifugation in processing murine and humanized mouse PBMCs and TILs for flow cytometry. Furthermore, we show how Laminar Wash improves flow cytometry data quality, in terms of debris removal and separation of immune cell subsets for both PBMCs and TILs.

**CONCLUSIONS:** Overall, these results show how Laminar Wash methodology assists in standardizing sample preparation for cytometric analysis, an important and unmet need in cancer immunotherapy discovery and manufacturing workflows.

## LAMINAR WASH SYSTEMS




## SEMI-AUTOMATED IMMUNOSTAINING WORKFLOW




**Figure 1.** Overview of CT26 syngeneic mouse model and Laminar Wash workflow. (a) CT26 colon tumor cells were transplanted subcutaneously to 6~8-week-old mice and established for approximately 2 to 3 weeks followed by the i.v. injection of PD-1 antibody or vehicle, respectively. Spleen and tumor samples were processed into single cell suspensions and analyzed for immune cell subsets by flow cytometry. (b) Overview of the sample preparation procedure using the Laminar Wash system. Dissociated tumor cells and splenocytes were transferred to a LW96 plate and washed on HT1000 during the staining procedure prior to flow cytometry.

## LAMINAR WASH IMPROVES RESOLUTION OF TUMOR INFILTRATING LYMPHOCYTES (TILs)




**Figure 2.** Flow cytometry gating strategy corresponding to centrifuge-processed (top) and Laminar Wash-processed (bottom) CT26 mouse tumor-naive samples. Samples were fixed prior to recording at an Attune NxT flow cytometry system. All percentages were gated from the parent population in the gating hierarchy.

## LAMINAR WASH IMPROVES RECOVERY OF VIABLE TILs



**Figure 3.** Splenocytes (a) and (b) and dissociated tumor cells (c) and (d) were washed with either centrifugation (manual 1 and 2) or the Laminar Wash system and viability measurements were compared. Statistical significance is reported among the manual methods vs Laminar Wash: ns = not significant, \*\* = P<0.01, \*\*\* = P<0.001, \*\*\*\* = P<0.0001. The values represent technical triplicates of the samples from an individual naïve and an individual challenged Balb/c mouse. Manually processed samples were handled by two different analysts.

## LAMINAR WASH REDUCES HANDS-ON TIME




**Figure 4.** Procedural workflows were considered either incubation (blue) or washing (cyan) steps and were subsequently timed accordingly. Laminar Wash has an additional initial settling time that is counted as a washing step in the conventional method. Incubation and settling are hands-free steps.

Mock (time-motion study)	Total Run Time, min:sec	Hands-On Time (% of total)	Walk-Away Time (% of total)
Centrifuge	25:06	10:06 (40.2%)	15:00 (59.8%)
HT2000	5:04	0:32 (10.5%)	4:32 (89.5%)

**Table 1.** A time-motion study comparing manual and Laminar Wash methodology was carried out, beginning after the antibody incubation step. Equivalent washes were performed: 3 runs of 5 min 1500 rpm spins on the centrifuge, 12 wash cycles with the Laminar Wash HT2000.


## LAMINAR WASH BASED CHARACTERIZATION OF VIRAL TRANSDUDED CAR T-CELLS



**Figure 5.** SupT1 cells infected with lentivirus were processed by iQue Screener Plus and read as 50µl per well. Laminar Wash shows improvements towards debris removal and viability staining.

**Figure 6.** Analysis for expression of lentivirus from SupT1 cells: Samples were stained and washed directly on Laminar Wash plates and the staining produced slightly higher fluorescence signal. Upon gating, the resolution for lentivirus expression was improved and well-to-well cleanliness was apparent. This performance indicates superior washing for SupT1 cells and consistency of readout. Further day-to-day performance was evaluated successfully.

## EFFICIENT WASHING OF IL-7 IN T-CELL STIMULATION KINETICS ASSAY



**Figure 7.** Strong TCR signals (caused by IL-7 mediated sensitization) cause developing thymocytes to activate apoptotic machinery (*negative selection*). Through the application of kinase inhibitors, those with targets that affect TCR signaling are able to override the process of negative selection: **T-cells are extremely sensitive to minute amounts of IL-7.** Centrifugation wash requires 3 rounds for complete dilution, plus stress on cells. Laminar Wash data demonstrates faster decay of pStat5 via efficient IL-7 removal.